

Lux Optical[®]

plum



GENERAL DESCRIPTION AND MATERIAL

- Eye wash solution
- 99.1% sterile water et 0,9 % sodium chlorure (corresponds to tears).
- Bottle of 500 ml

SPECIFICS ADVANTAGES

- Ergonomic design for easy eye washing
- Anti-dust
- Clear use instructions
- Available alone or in care station

MAIN USED FIELDS

- For classic eye washing against dust
- User information through pictograms and mirror in the care station
- Shelf life : 3 years after manufacturing date mentioned on the bottle
- For manufacturing companies, laboratories, ...
- Do not place it in dusty or dirty environments
- Open bottles or remaining liquids must not be reused.

PACKAGING

- Carton of 12 bottles
- Store between 5 and 35 °C

STANDARDIZATION

All PLUM eye wash items are **CE marked** according to regulations on the current **medical devices** (European Directive 93/42/EEC, class I). They are sterile so that the eye is not infected by foreign germs during the rinsing.

WARNING !

Whatever the criticality of the accident, the pressure on the bottle during the washing must be soft and continu.

In case of ocular damages, consult a doctor. The washing must be performed until the come of the doctor.